

HO-PEG-NH-Fmoc (MW 1000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	HO-PEG-NH-Fmoc (MW 1000) is a PEG derivative featuring an Fmoc protecting group and a hydroxyl (OH) reactive group. The Fmoc group is a commonly utilized amino protecting group that safeguards the amino group from unintended reactions and can be removed when a specific coupling reaction is desired.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481