

BOC-NH-PEG-FA (MW 10000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	BOC-NH-PEG-FA (MW 10000) [BOC-NH-PEG-Folate (MW 10000)] is a PEG derivative modified with Folic acid. Folic acid (-FA) has a high affinity for folate receptors, which can be utilized for targeted drug delivery to cell membrane receptors.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481