

XL-3156

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	XL-3156 is a potent, selective, and cross-species inhibitor of cGAS. It simultaneously occupies allosteric and orthosteric sites, stabilizing the closed conformation of the activation loop to inhibit the interaction and phase separation of cGAS with DNA. XL-3156 is applicable for research into autoimmune diseases and inflammation.
Targets(IC50)	cGAS

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481