

FRETS-VWF73 TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from moisture**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	FRETS-VWF73 TFA is a 73-amino acid peptide that serves as a fluorescent substrate for ADAMTS13 detection (Ex=340 nm; Em=450 nm).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481