

18:0 SquareFluor(629/648) PE ammonium

Chemical Properties

CAS No. : 3077426-99-1

Formula:

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	18:0 SquareFluor(629/648) PE ammonium is a fluorescent lipid used for visualizing cell membranes (Ex/Em = 629/648 nm).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481