

Targeting the bacterial sliding clamp peptide 46 TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Targeting the bacterial sliding clamp peptide 46 TFA is a short peptide that specifically targets the bacterial sliding clamp (SC) and inhibits SC-dependent DNA synthesis.
Targets(IC50)	Antibacterial

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481