

PcrA spirase

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PcrA spirase is an ATP-dependent helicase that unwinds double-stranded DNA (dsDNA) from the 3' to 5' direction. It binds to dsDNA and separates it into two single-stranded DNA (ssDNA) molecules.
Targets(IC50)	DNA/RNA Synthesis

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481