

RPX 7546

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	RPX 7546 is a potent inhibitor of metal- $\beta$ -lactamases (MBLs), demonstrating IC50 values of 1.4 $\mu$ M for IMP-1 and 48 $\mu$ M for IMP-78. Additionally, RPX 7546 exhibits antibacterial activity.
Targets(IC50)	Antibacterial

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481