

1,2-Dioleoyl-sn-glycero-3-PA ammonium salt

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	1,2-Dioleoyl-sn-glycero-3-PA (ammonium salt) is a phospholipid characterized by the incorporation of the long-chain fatty acid oleic acid at the sn-1 and sn-2 positions. It is utilized in the synthesis of artificial membranes such as liposomes.
Targets(IC50)	Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481