

RS-79948-197 hemihydrate

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	RS-79948-197 hemihydrate is a non-imidazoline α 2-adrenoceptor antagonist. It exhibits Kd values of 0.42 nM, 0.18 nM, and 0.19 nM for rat α 2A, α 2B, and α 2C receptors, respectively, and 0.60 nM, 0.46 nM, and 0.77 nM for human α 2A, α 2B, and α 2C receptors.
Targets(IC50)	Adrenergic Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481