

## Boc-Gly-Lys-Arg-AMC diacetate

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Boc-Gly-Lys-Arg-AMC diacetate is a polypeptide that can be identified through polypeptide screening. This screening method is a research tool primarily using immunoassays to gather active polypeptides. It is useful in studying protein interactions, conducting functional analysis, screening antigen epitopes, and is particularly significant in the field of active molecule research and development.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481