

TE-8105

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

|               |   |
|---------------|---|
| Description   | TE-8105 is a GLP-1 receptor agonist that demonstrates prolonged and enhanced efficacy in models of diabetes, obesity, and non-alcoholic steatohepatitis (NASH). |
| Targets(IC50) | Glucagon Receptor   |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481