

c-Met ligand-Linker Conjugate 1

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	c-Met ligand-Linker Conjugate 1 is a Target Protein Ligand-Linker Conjugate that consists of a c-Met ligand and a PROTAC linker, capable of recruiting E3 ligase. This compound is used in the synthesis of PROTACc-Met degrader-5.
Targets(IC50)	c-Met/HGFR,Target Protein Ligand-Linker Conjugates

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481