

## K-Ras ligand-Linker Conjugate 6 formate

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	K-Rasligand-Linker Conjugate 6 formate includes a target protein K-Ras ligand (K-Ras recruitment motif) and a PROTAC linker, capable of recruiting E3 ligases such as VHL, CRBN, MDM2, and IAP. It is utilized in synthesizing PROTACK-Ras Degrader-1, a potent PROTACK-Ras degrader, achieving a K-Ras degradation rate of $\geq 70\%$ in SW1573.
Targets(IC50)	Target Protein Ligand-Linker Conjugates

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481