

DSPE PEG2k-OH ammonium

Chemical Properties

CAS No. : 3056014-08-2

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE PEG2k-OH ammonium is a hydroxyl-terminated PEG phospholipid. Its hydrophobic tail can encapsulate and aggregate other hydrophobic drugs, while the hydroxyl end is suitable for further reactions. DSPE PEG2k-OH ammonium is applicable for the preparation of liposomes or lipid nanoparticles.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481