

PI4K-IN-2

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PI4K-IN-2 (Compound 30) is a highly selective inhibitor of type III phosphatidylinositol 4-kinase beta (PI4KB) with an IC50 value of 0.015 μ M. It holds potential for research in human rhinovirus (HRV).
Targets(IC50)	PI4K, Virus Protease

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481