

(S,R,S)-AHPC-PEG3-ADIBO-adipic acid

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	(S,R,S)-AHPC-PEG3-ADIBO-adipic acid (Compound 2a) is a conjugate of the PROTAC ligand (Compound 1d) and linker. This compound can also couple with an aptamer to construct a PROTAC.
Targets(IC50)	Others,PROTACs

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481