

## Canine KP-10

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from moisture**  
**Powder: -20°C for 3 years | In solvent: -80°C for 1 year**  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Canine KP-10 is a potent kisspeptin that induces a significant gonadotropin and estradiol response in female dogs during estrus.
Targets(IC50)	Kisspeptin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481