

HPK1 ligand 4

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	HPK1 ligand 4 is a PROTAC (proteolysis-targeting chimeras) target protein ligand used in the synthesis of PROTACs, such as PROTACHPK1 Degradator-5. PROTACHPK1 Degradator-5 is a potent and orally active HPK1 PROTAC degrader demonstrating antitumor activity.
Targets(IC50)	MAPK,Ligands for Target Protein for PROTAC

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481