

## Anti-MRSA agent 34

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-MRSA agent 34 (Compound 6) is an antibacterial agent that targets methicillin-resistant Staphylococcus aureus (MIC=14 µg/mL). It holds potential for research into microbial infections.
Targets(IC50)	Antibacterial

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481