

## BRD4 ligand 8

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Keep away from direct sunlight

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	BRD4 ligand 8 (Compound J558) is an inhibitor of BRD4. This compound can be utilized in the synthesis of PROTACBRD4 Degradar-33.
Targets(IC50)	Epigenetic Reader Domain,Ligands for Target Protein for PROTAC

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481