

[2(R,S)-2-Sulfanylheptanoyl]-Phe-Ala-CoA

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	[2(R,S)-2-Sulfanylheptanoyl]-Phe-Ala-CoA ([2(R,S)-2-Sulfanylheptanoyl]-Phe-Ala-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481