

## Sp-dGTPαS tetrasodium

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Sp-dGTPαS tetrasodium is a sulfur-containing isomer of dGTP. It acts as an inhibitor of SAMHD1 triphosphate hydrolase activity, with a $K_i$ value of 820 nM.
Targets(IC50)	Nucleoside Antimetabolite/Analog

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481