

## Hepoxilin A3 methyl ester

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Hepoxilin A3 methyl ester is the ester form of HxA3. It effectively induces the release of intracellular calcium ions, subsequently facilitating calcium influx in suspended human neutrophils.
Targets(IC50)	Calcium Channel

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481