

DSPE PEG(2000)-N-Cyanine 7

Chemical Properties

CAS No. : 2315262-01-0

Formula:

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE PEG(2000)-N-Cyanine 7 is a phospholipid characterized by excellent biocompatibility and significant amphiphilicity, making it a primary and suitable form or excipient in pharmaceutical formulations to achieve enhanced therapeutic outcomes. DSPE PEG(2000)-N-Cyanine 7 can be used in drug delivery research.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481