

## Pomalidomide-PEG3-Cys

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Pomalidomide-PEG3-Cys is an E3 ubiquitin ligase ligand-linker conjugate, consisting of a CRBN ligand Pomalidomide and a three-unit PEG linker. It can be utilized in the synthesis of PROTAC BRD4 Degrader-33.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481