

## Maleimide-Ph(3,5-F)-PEG4-Val-Ala-MF-6

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Maleimide-Ph(3,5-F)-PEG4-Val-Ala-MF-6 is a drug-linker conjugate for ADCs. It consists of a topoisomerase 1 inhibitor (Topi MF-6) and a cleavable linker (Maleimide-Ph(3,5-F)-PEG4-Val-Ala). This compound is utilized in the synthesis of antibody-drug conjugates (ADCs).
Targets(IC50)	Drug-Linker Conjugates for ADC,Topoisomerase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481