

## LAG-3 Hit II

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	LAG-3Hit II is an inhibitor of LAG-3 with a KD of 4.32 $\mu$ M. It is applicable in research related to cancers such as ovarian cancer, colon adenocarcinoma, and melanoma.
Targets(IC50)	LAG-3

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481