

SRC-1 NR box peptide acetate

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	SRC-1 NR box peptide acetate is a bioactive peptide, consisting of a 14-amino acid fragment derived from the steroid receptor coactivator SRC-1 NR II. This compound is useful for investigating the regulatory mechanisms of estrogen receptor ligands.
Targets(IC50)	Estrogen Receptor/ERR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481