

AAK1-IN-11

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	AAK1-IN-11 is an inhibitor of AAK1 (adaptor associated kinase 1) with an IC50 of 2 nM, and it exhibits antiviral activity.
Targets(IC50)	AAK1 (AP2 associated kinase 1)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481