

AH081

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	AH081 (Compound 38) is a CK1 $\delta/\epsilon$ PROTAC degrader and serves as the negative control for AH078. It maintains inhibitory activity against CK1 $\delta/\epsilon$ but lacks degradation activity due to the use of an inactive stereoisomer of the VHL ligand.
Targets(IC50)	Casein Kinase,PROTACs

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481