

DSPE-PEG2k-tri Mannose ammonium

Chemical Properties

CAS No. : 3056016-73-7

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG2k-tri Mannose ammonium is a tri-mannose-DSPE-conjugated PEG derivative used in the synthesis of lipid nanoparticles for drug delivery.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481