

G-Pen-GRGDSPCA TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	G-Pen-GRGDSPCA TFA is an $\alpha\beta3$ inhibitory RGD peptide. It inhibits neointimal hyperplasia and luminal narrowing following vascular injury and regulates the migration and proliferation of smooth muscle cells. This compound can be utilized in studies related to mechanisms of vascular injury repair.
Targets(IC50)	Integrin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481