

## XIAP ligand 5

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	XIAPligand 5 is an XIAP antagonist and serves as an E3 ligase ligand (E3 ligase Ligand). This compound is utilized in the synthesis of PROTACGNE-1567.
Targets(IC50)	Ligands for E3 Ligase,IAP

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481