

Rp-TTPαS tetrasodium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Rp-TTPαS (tetrasodium) is the Rp isomer of TTPαS (tetrasodium). It can bind to HIV-1 reverse transcriptase (HIV-1RT) and exhibits a significant thio-phosphate element effect.
Targets(IC50)	HIV Protease

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481