

SARS-CoV-2-IN-118

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	SARS-CoV-2-IN-118 (compound A.13) is an inhibitor that targets the translation initiation element SL1 of SARS-CoV-2 RNA.
Targets(IC50)	SARS-CoV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481