

## Chol-PEG2000 fluorescein

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Chol-PEG2000 fluorescein is a fluorescent lipid composed of a cholesterol-PEG2000 unit conjugated with fluorescein. This compound is suitable for preparing lipid nanoparticles for drug delivery research [Ex/Em = 480/520 nm].
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481