

DOPE-PEG(2000) Carboxylic acid sodium

Chemical Properties

CAS No. : 2342575-87-3

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DOPE-PEG(2000) Carboxylic acid sodium is a PEG polymer containing unsaturated DOPE phospholipids. It is utilized in the preparation of liposomes or other nanoparticles for drug delivery applications.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481