

MG-HuR2

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	MG-HuR2 is a molecular glue degrader that targets the oncogenic RNA-binding protein HuR, with an IC50 of 0.5 μ M. It holds potential for research on HuR-overexpressing malignancies, such as breast cancer.
Targets(IC50)	Molecular Glues,HuR,PROTACs

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481