

RIPK1-ligand-2

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	RIPK1-ligand-2 is an inhibitor of RIPK1, and it can be used in the synthesis of PROTAC [RIPK1]Degradar-1.
Targets(IC50)	Ligands for Target Protein for PROTAC,RIP kinase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481