

CES

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	CES (CESPLLSEC) is a homing peptide that targets synaptosome-associated protein 25 (SNAP25). It specifically homes to intracranial U87MG and WT-glioblastoma (GBM) models. CES can be utilized as a targeting peptide for drug delivery in cancer research.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481