

## DOPE-PEG(2000)-HALO-TAG ammonium

### Chemical Properties

CAS No. : 2342575-90-8

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	DOPE-PEG(2000)-HALO-TAG ammonium is a PEG polymer incorporating unsaturated DOPE phospholipids. It is suitable for the preparation of liposomes or other nanoparticles employed in drug delivery.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481