

IL-23 cyclic peptide inhibitor 281

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	IL-23cyclic peptide inhibitor 281 is a cyclic peptide that acts as an IL-23 inhibitor. It is applicable for research in inflammatory diseases, such as inflammatory bowel disease (IBD).
Targets(IC50)	Interleukin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481