

2'-5' Oligoadenylate Synthetase 1, Human

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	2'-5' Oligoadenylate Synthetase 1, Human is an interferon-induced enzyme that is activated in the presence of double-stranded RNA, promoting the oligomerization of ATP into 2',5'-linked oligoadenylates (2-5A). It plays a crucial role in the inflammatory immune response.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481