

RIPK1 ligand-Linker Conjugate-1

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	RIPK1 ligand-Linker Conjugate-1 is a Target Protein Ligand-Linker Conjugate that consists of a RIPK1 ligand and a PROTAC linker, designed to recruit E3 ligase. It is utilized in the synthesis of PROTACRIPK1 Degradator-1.
Targets(IC50)	RIP kinase, Target Protein Ligand-Linker Conjugates

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481