

CDK12-Cyclin K ligand-1

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	CDK12-Cyclin K ligand-1 is a ligand for PROTAC degraders, designed to bind with E3 ligase. It is utilized in synthesizing PROTAC degraders, such as [PP-C8].
Targets(IC50)	Ligands for E3 Ligase,Ligands for Target Protein for PROTAC

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481