

DHODH-IN-30

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DHODH-IN-30 (compound 4g) is a potent inhibitor of human dihydroorotate dehydrogenase (hDHODH), with an IC50 of 32 nM. Additionally, DHODH-IN-30 offers protective effects against oxidative stress.
Targets(IC50)	Dehydrogenase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481