

RB-PEG-SH (MW 400)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	RB-PEG-SH (MW 400) (Rhodamine B-PEG-SH (MW 400)) is a fluorescent dye composed of Rhodamine B, PEG, and a thiol group (SH). The SH group is highly reactive, capable of forming stable thioether bonds with various functional groups, such as maleimide (Ex/Em = 546/610 nm).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481