

Arthrofactin TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Keep away from moisture

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Arthrofactin (TFA) is an effective lipopeptide biosurfactant from Pseudomonas sp. MIS38. Its production is associated with various ATP-dependent active transport systems.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481