

Glutamate Dehydrogenase, Bovine Liver

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Glutamate Dehydrogenase, Bovine Liver (EC 1.4.1.4), catalyzes the reversible oxidative deamination of glutamate, producing α -ketoglutarate and ammonia.
Targets(IC50)	Dehydrogenase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481